BÁSCULA DE BAÑO ELECTRÓNICA – MANUAL DE INSTRUCCIONES
ELECTRONIC BATHROOM SCALE-USER MANUAL
BALANCE DE SALLE ÉLECTRONIQUE – MANUEL D’INSTRUCTIONS
ESCALA DE BANHEIRO ELECTRÓNICA- MANUAL DE INSTRUÇÕES

PB 2232

Sonifer, S.A.
Avenida de Santiago, 86
30007 Murcia
España
E-mail: sonifer@sonifer.es / web: www.orbegozo.com
MADE IN P.R.C.

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ENGLISH**SAFETY WARNINGS**

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing.

1. This appliance can be used by children aged from 8 years and above and persons with reduced physical, sensory or mental capabilities or lack of experience and knowledge if they have been given supervision or instruction concerning use of the appliance in a safe way and understand the hazards involved. Children must never play with the appliance. Cleaning and user maintenance must never be carried out by children without supervision.
2. Children should be supervised to ensure that they do not play with the appliance.

3. Keep the appliance out of the reach of children less than 8 years old.
4. **WARNING:** In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.
5. Do not use the unit if it is not working properly.
6. Do not handle the appliance with wet hands.
7. Never immerse the appliance in water or any other liquid.
8. Make sure the appliance has been unplugged before cleaning.
9. This appliance must be installed following the national regulations for electrical installations.
10. This appliance is for household use only.

1. During Measuring

- A. The device is designed to allow auto step on. Please always ignore the first reading and only start taking effective record from the second time weighing.
- B. Always remove your and , wipe clean your feet before using.
- C. Step on the weighing platform gently.

- D. It is advisable to take measurement at the same time of the day.
- E. Measured result could be misleading after intensive exercise, excessive dieting or under extreme dehydration condition.
- F. Always weigh/use the scale on a hard and flat surface. When you are measuring, please keep still.
- G. The measured data of the following people may have deviation:
Children under 10 or adult over 100 (can only use the device under normal weighing mode)
Adult over 70
Body building athletes or other occupational athletes.
- H. People with the following conditions should not use the device:
Pregnant women
Those with symptom of edema
Those on dialysis treatment
Those who use heart pacemaker or with other implanted medical device.

2. Safety Warning

- A. The Fat% Hydration% content evaluated by the Fat Monitor Scale is only for reference (not for medical purpose). If your Fat% Hydration% content is over/below normal level, please consult your doctor for more advice.
- B. Keep away from strong electro-magnetic field when using the scale.
- C. It could be slippery if you weigh on a wet surface. Moreover, the weighing platform surface is glazed, never step on the edge of the scale and do keep balance during your weighing.
- D. **WARNING:** In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

Weighing Operation


1. Insert the batteries 2 x 1.5V AAA into the battery compartment with correct polarity.
2. Turn on scale, choose a measuring unit with up/down buttons
3. Place the scales on a firm, flat and even surface to ensure accuracy.
4. Step on the scales platform to switch on the product, stand still while the scales compute your weight.
5. Your weight will be displayed and will flash twice before it is fixed.
6. When you step off the scales it will automatically switch off after 15 seconds.
7. If no weighing is carried out while 0.0 is displayed the scales will automatically switch off after 15 seconds.

Remarks: This initialization process as below must be repeated if the scale is moved. At all other times step straight on the scale.

- Step 1. Press the platform centre and remove your foot.
- Step 2. "0.0" will be displayed
- Step 3. The scale will switch off and now it is ready for use.

This initialisation process must be repeated if the scale is moved. At all other times step straight on the scale.

- Step 4. Gently step on the scale. Put your feet evenly on the scale without moving and wait whilst the scale compute weight. Your correct weight will then be displayed in the following way: 80.3 Kg.

Body Fat Measurement:

In order to ensure accurate results users should enter the required personal parameters:

age, sex and height.

These scales have the facility to store the personal parameters of 10 people.

The technology used in these scales, Bio-electrical impedance Analysis (BIA) is more technologically advanced than the older Body Mass Index (BMI) method where fat % is calculated using only height and weight. BIA technology uses 5 parameters to get a more accurate resolution for the analysis (weight and biological resistance which are automatically measured by the scales together with the three personal parameters entered manually).

Mode selection and parameter setting:

Step on the scale platform to switch on

Press UP or DOWN to switch between Body Parameter 0-9. Press SET to confirm.

Press UP or DOWN to select male or female. Press SET to confirm

Press UP or DOWN to indicate height Press SET to confirm

Press UP or DOWN to indicate age. Press SET to confirm

The values are flashed and recorded.

Step on the scales after 0.0 appears, First the weight is measured and recorded, then is displayed, indicating that the scales are analysing, please do not alight. Once finished, the weight and body fat percentage, water percentage, bone weight, muscle percentage and KCAL and BMI will be displayed alternately for 2 times before the scales automatically turn off. Refer to the body fat reading against the (Body Fat & Water Ratio Chart) to find out the body status. please find the screen display order as following

Weight display

Body fat percentage display

Body water percentage display


Body muscle percentage

Body bone weight display

KCAL display

Body fat ratio Chart(%)

AGE	FEMALE(%)				MALE(%)			
	underweight	healthy	overweight	obese	underweight	healthy	overweight	obese
6-39	<21	21~33	33~39	>39	<8	8~19	19~25	>25
40-59	<23	23~35	35~40	>40	<11	11~22	22~28	>28
60-100	<24	24~36	36~42	>42	<13	13~25	25~30	>30

Due to the natural difference between people this chart should be used for reference purposes only.

Water Ratio Chart & Muscle Ratio Chart (%)

AGE	FEMALE (%)				MALE (%)			
	Low	Healthy	High	Body Muscle	Low	Healthy	High	Body Muscle
6-15	<57.0	57.0-67.0	>67.0	>34	<58.0	58.0-72.0	>72.0	>40
16-30	<47.0	47.0-57.0	>57.0		<53.0	53.0-67.0	>67.0	
31-60	<42.0	42.0-52.0	>52.0		<47.0	47.0-61.0	>61.0	
61-100	<37.0	37.0-47.0	>47.0		<42.0	42.0-56.0	>56.0	

Bone Weight Chart (Kg)

	FEMALE			MALE		
	<45 Kg	45-60 Kg	>60 Kg	<60 Kg	60-75 Kg	>75 Kg
Ref.bone weight	-2.6 Kg	-3.2Kg	-8Kg	-3.4 Kg	-3.8 Kg	-8Kg

KCAL Chart

AGE	FEMALE		MALE	
	WEIGHT(KG)	KCAL	WEIGHT(KG)	KCAL
6-17	50	1265	60	1620
18-29	55	1298	65	1560
30-49	60	1416	70	1680
50-69	60	1302	70	1561
70-100	50	1035	65	1397

Trouble shooting

1. When the LCD screen displays "LO", you should install a new battery.
2. When the LCD displays "ERR", it indicates that the maximum weight has been exceeded.
3. When the LCD screen displays Err 2 ,it shows an error in the body fat percentage

Use and care

- Ensure the sales are on a flat, even and firm surface .Avoid carpet or soft surfaces.
- Always weight yourself with the same scales placed on the same surface .Uneven floors may affect the reading.The surface of the scales will be slippery if wet. Keep it dry!

- Weight yourself without clothing and footwear, before meals and always at the same time of day.
- For better accuracy, it is not advised to take the measurement within 2 hours of waking up.
- Keep the scales in a dry place to protect the electronic components.
- Clean after use with a lightly damp cloth, do not use solvents or immerse the product in water.
- Avoid overloading the scales; otherwise, permanent damage may occur.
- Do not disassemble the scales other than to replace the battery.
- Do not store the scales in an upright position when not in use as this will drain the battery.
- Avoid excessive impact or vibration to the unit, such as dropping the scales onto the floor, or objects onto the scales.

Disposal of old electrical appliances.

The European directive 2012/19/UE on Waste Electrical and Electronic Equipment (WEEE), requires that old household electrical appliances must not be disposed of in the normal unsorted municipal waste stream. Old appliances must be collected separately in order to optimize the recovery and recycling of the materials they contain, and reduce the impact on human health and the environment.

The crossed out "wheeled bin" symbol on the product reminds you of your obligation, that when you dispose of the appliance, it must be separately collected

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CERTIFICATE OF GUARANTEE

We hereby certify that the equipment stated in this Guarantee Card has gone through strict quality control checks and bears a full Guarantee of 24 months, which shall come in force from the date of purchase.

During the 24-month guarantee period, the present certificate covers repairs, workmanship and spare-parts for all manufacturing defects, however this guarantee does not cover those defects that have been caused due to the following: defects of installation, blows, bad use, breakage of plastic parts and external causes not related to the equipment. This guarantee does not include installation and maintenance.

This certificate of guarantee is cancelled on those accessories and components which are worn due to continuous use. Those which are perishable, such as plastic components, glass, lamps, filters, rubber or enamelworks which get worn due to wrongful use or reaction to outside agents such as heating, wáter or chemical products.

The present certificate of guarantee shall be nullified if the equipment is handled by any unauthorised personnel.

Customer's rights:

The seller should be responsible over the customer in the event of lack of conformity when the product is purchased. It is recognised the right of the customer to get the product repaired, substituted, decreased its price or the resolution of the contract.

The present certificate is only valid if it is accompanied by the receipt or invoice of purchase.

MODEL _____ Nº _____
PURCHASER'S NAME _____
ADDRESS _____ Nº _____
TOWN _____ ZIP CODE _____
DATE OF PURCHASE _____
RETAILER'S STAMP _____

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